


SURG-I-LOOP® Vessel Loops



Designed to assist surgeons in complicated surgical fields by providing retraction, occlusion and identification of arteries, veins, nerves and ureters.

- **Unique design** – provides access to complicated surgical fields
- **Versatile** – can retract and occlude arteries, veins, nerves and ureters
- **Durable** – retains preset tension
- **Radiopaque**
- **Not made with natural rubber latex**

Cat. No.	Description
1001-75	Red, Mini
1001-76	Red, Maxi
1001-77	Blue, Mini
1001-78	Blue, Maxi
1001-78A	Blue, Super Maxi
1001-79	Yellow, Mini
1001-80	Yellow, Maxi
1001-81	White, Mini
1001-81M	White, Micro
1001-82	White, Maxi

Length 16.5" (42 cm)

2 per sterile package, 10 packages per box.

	Super Maxi (5 x 1.27 mm)
	Maxi (2.5 x 1.02 mm)
	Mini (1.5 x .88 mm)
	Micro (0.5 x .26 mm)

SCANLAN® SURG-I-LOOP® PLUS

SURG-I-LOOP® Plus features a solid stainless steel tip for fast, accurate placement / passing in complicated surgical fields. Ideal for minimally-invasive cardiac surgery. Radiopaque. Not made with natural rubber latex.

Cat. No.	Description
1001-81P	Blunt needle (1.2 mm) (light green)
1001-82P	Semi-sharp needle (1.2 mm) (light orange)

2 per sterile package, 10 packages per box.

