




VATS and Open Thoracotomy Procedures

- May be useful for:
 - Chest wall resections
 - First rib resections
 - Resecting the costal cartilage during some Pectus procedures
- Smooth feeling cutting mechanism provides clean, sharp cuts
- Safety - the blade is centered on the footplate minimizing the possibility of inadvertently cutting tissue surrounding the rib or costal cartilage
- 0-6 cm ruler on the distal end of the instrument




9909-900



9909-902



9909-904
9909-905



Cat. No.	Blade Height	Cutting Width
9909-900	10 mm	16 mm
9909-902	15 mm	25 mm
9909-904	20 mm	25 mm

SCANLAN®
Dennis Rib Cutter

Cat. No.	Description	Shaft	Working Length	Total Length
9909-900	10 mm blade 16 mm cutting width Straight shaft	10 mm	11 1/2" (29 cm)	13 1/2" (34 cm)
9909-902	15 mm blade 25 mm cutting width Straight shaft			
9909-904	20 mm blade 25 mm cutting width Straight shaft			
9909-905	20 mm blade 25 mm cutting width Curved shaft			

Please Note...
Each Scanlan instrument is individually handcrafted to provide the surgeon with the highest quality, longest-lasting, precision instrument possible. Reflecting this unique process, all stated dimensions are approximate.

Please check with your Scanlan Representative for Firm Sale status. Custom modifications available upon request.
For a complete listing of models visit... www.scanlaninternational.com

