

CHROMAPRO™ LIGHT SOURCE

CHROMAPRO LED LIGHT SOURCE



SEE BETTER  
SEE FARTHER  
SEE WHY



SINCE  
1971

CELEBRATING  
FIVE DECADES  
OF ILLUMINATION

## CHROMAPRO™ LED LIGHT SOURCE

### LED LIGHT SOURCE

- 50-Watt LED Illumination
- Illumination Range with AtoN Fiber Optic Headlight -70,000 LUX at 16" (40 cm)
- LED Expected Lifespan > 30,000 hours
- Analog Dimmer Controls Output from 0% to 100%
- Universal Power Supply
- Optimizes Light Transmission through Cable
- Small and Powerful

### CHROMAPRO PRODUCT CODES:

BFW 5500	ChromaPRO LED Light Source
BFW 1200	AtoN Headlight with Fiji Headband and 10' Cobra Coil Fiber Optic Cable
BFW 1255	AtoN Headlight with Fiji Headband and 10' Standard Fiber Optic Cable
BFW 1260	AtoN Headlight with Fiji Headband and 7' Standard Fiber Optic Cable

### EXCELLENT FOR USE IN THE FOLLOWING DISCIPLINES

ARTHROSCOPY  
COSMETIC SURGERY  
ENT  
EXAM  
ORAL/MAXILLOFACIAL  
OUTPATIENT SURGERY  
PLASTIC SURGERY  
PRIVATE PRACTICE



[BFWinc.com](http://BFWinc.com)

800-717-4673 | +1-502-899-1808

1550 SAMCO ROAD

RAPID CITY, SD, 57702 USA

**ENGINEERED AND ASSEMBLED IN THE U.S.A.** © 2024 BFW, Inc.  
ALL TECHNICAL SPECIFICATIONS ARE SUBJECT TO CHANGE WITHOUT NOTICE  
**CHROMAPRO TESTED AND CERTIFIED TO:** IEC/EN 60601-1:2006, UL 60601-1:2003, CS C22.2 No. 601.1 1990 (R2005), and IEC 60601-1-2:2007