

DOVER™ | SURGICAL & MEDICAL

PORTABLE LED HEADLIGHT SYSTEM



SEE BETTER
SEE FARTHER
SEE WHY




SINCE
1971

CELEBRATING
FIVE DECADES
OF ILLUMINATION

DOVER™ EXAM: MEDIUM INTENSITY

PORTABLE LED HEADLIGHT WITH FIXED SPOT, POWERED BY A 3-CELL BATTERY

- Illumination Range: ~30,000 LUX at 16" (40 cm) to ~105,000 LUX at 10" (25 cm)
- Fixed Spot Diameter: 3" (7.6 cm) at a working distance of 16" (40 cm)
- Color Temperature: 7,000 K
- Intensity Settings: 5 Levels
- 3-Cell Battery
- Battery Duration at Full Intensity: 10 hours
- Battery Expected Life Span: > 300 Charges
- LED Expected Life Span: > 10,000 hours
- Optic Weight: 0.45 oz (12.7 g)
- Battery Weight: 6.5 oz (184 g)



DOVER PORTABLE LED HEADLIGHT SYSTEM INCLUDES:

- Dover High-Intensity Headlight
- Choice of: Fiji or Key West Headband, or Optic Clip
- Choice of One or Two 3-Cell Battery
- Charging Cradle, Power Supply
- Holster, Linkage Kit

DOVER PORTABLE LED HEADLIGHT SYSTEM PRODUCT CODES:

BFW 9030	Dover Headlight System with Fiji Headband with One 3-Cell Battery
BFW 9030B	Dover Headlight System with Fiji Headband and Two 3-Cell Batteries
BFW 9032	Dover Headlight System with Optic Clip with One 3-Cell Battery
BFW 9032B	Dover Headlight System with Optic Clip and Two 3-Cell Batteries
BFW 9034	Dover Headlight System with Key West Headband with One 3-Cell Battery
BFW 9034B	Dover Headlight System with Key West and Two 3-Cell Batteries
BFW 1552C	Dover Optic Clip Adapter for Face Shields and Loupes
BFW 5440	Dover Protective Carrying Case, molded interior provides space for all components
BFW 9022ACC	Dover Charger
BFW 9035	Dover 3-Cell Battery
BFW 9036	Dover Holster

EXCELLENT FOR USE IN THE FOLLOWING DEPARTMENTS

HOSPITAL	ENT	SURGERY CENTER	ORTHOPEDIC SURGERY	
	GENERAL SURGERY		PLASTIC SURGERY	
	CARDIOSURGERY			
	NEUROSURGERY			
	ORTHOPEDIC SURGERY		OFFICE	DENTAL
	PEDIATRIC SURGERY			EXAM
SPINE SURGERY	ORAL/MAXILLOFACIAL			
THORACIC SURGERY				
TRANSPLANT SURGERY				
VASCULAR SURGERY				



BFWinc.com

800-717-4673 | +1-502-899-1808

1550 SAMCO ROAD

RAPID CITY, SD, 57702 USA

ENGINEERED AND ASSEMBLED IN THE U.S.A. © 2024 BFW, Inc.

ALL TECHNICAL SPECIFICATIONS ARE SUBJECT TO CHANGE WITHOUT NOTICE

DOVER TESTED AND CERTIFIED TO: IEC 60601-1 Edition 3, IEC 60601-1-2 Edition 4, ANSI/AAMI ES60601-1 (2005 + C1:09 + A2:10), CAN/CSA C22.2 NO 60601-1 Edition 2 (2008), CE