

NEW!

MAUI CORDLESS HEADBAND WITH BRISTOL-PLUS OPTIC



**SEE BETTER
SEE FARTHER
SEE WHY**




**SINCE
1971**

**CELEBRATING
FIVE DECADES
OF ILLUMINATION**

NEW! MAUI BRISTOL-PLUS™

CORDLESS LED SURGICAL HEADLIGHT WITH ADJUSTABLE SPOT, POWERED BY A HEADBAND MOUNTED BATTERY

- Illumination Range: ~40,000 LUX at 16" (40 cm) to ~132,000 LUX at 10" (25 cm)
- Adjustable Spot Diameter: 1"-3.5" (2.54-8.89 cm) at a working distance of 16" (40 cm)
- Color Temperature: 7,000 K
- Intensity Controls Output from 0% to 100%
- Battery Expected Life Span: > 300 Charges
- LED Expected Life Span: > 50,000 hours
- Battery Run Time: 3 hours at full intensity, +1 hour of decreasing intensity
- Battery Charge Time: 3 hours
- Optic Weight: 1.2 oz (34.02 g)
- Battery Weight: 1.8 oz (51 g)
- Total Weight: 9.9 oz (280.66 g)



MAUI BRISTOL-PLUS CORDLESS LED HEADLIGHT SYSTEM INCLUDES:

- Bristol-Plus Optic on Maui Headband
- Two Mini Batteries
- Dual-Bay Charging Cradle, Power Supply
- Linkage Kit

MAUI BRISTOL-PLUS CORDLESS LED HEADLIGHT PRODUCT CODES:

- BFW 9040PM** Maui Bristol-Plus Cordless Headlight System with Two Mini Batteries and Dual-Bay Charger
- BFW 9035M** Maui Bristol-Plus Mini Battery
- BFW 9023M** Dual-Bay Charger with Wall Power Supply for USA, EU, UK, Canada, Mexico, Australia
- BFW 9023MDPS** Dual-Bay Desktop Charger with Power Supply for all other countries
- BFW 5445** Protective Carrying Case

EXCELLENT FOR USE IN THE FOLLOWING DISCIPLINES

HOSPITAL

- ENT
- GENERAL SURGERY
- CARDIOSURGERY
- NEUROSURGERY
- ORTHOPEDIC SURGERY
- PEDIATRIC SURGERY
- SPINE SURGERY
- THORACIC SURGERY
- TRANSPLANT SURGERY
- VASCULAR SURGERY

SURGERY CENTER

- ORTHOPEDIC SURGERY
- PLASTIC SURGERY

OFFICE

- DENTAL
- EXAM
- ORAL/MAXILLOFACIAL



BFWinc.com

800-717-4673 | +1-502-899-1808

1550 SAMCO ROAD

RAPID CITY, SD, 57702 USA

ENGINEERED AND ASSEMBLED IN THE U.S.A. © 2024 BFW, Inc.

ALL TECHNICAL SPECIFICATIONS ARE SUBJECT TO CHANGE WITHOUT NOTICE

MAUI BRISTOL-PLUS TESTED AND CERTIFIED TO: IEC 60601-1 EDITION 3, IEC 60601-1-2 EDITION 4, AAMI ES60601-1 EDITION 1, CSA C22.2 NO 60601-1:08 EDITION 2, CE