


SINGLE-USE



SCANLAN® SUTURE-STATT®

Suture Tags

Designed to facilitate the grasping and tagging of suture during ophthalmic or delicate procedures.

- **Secure** – firmly holds the finest suture without damage
- **Keeps sutures out of the operative field**
- **Lightweight** – only 0.5 grams
- **Radiopaque**
- **Not made with natural rubber latex**

| Cat. No. | Description |
|----------|-------------|
|----------|-------------|

| | |
|----------|-------------|
| 1001-520 | Suture Tags |
|----------|-------------|

2 per sterile package, 5 packages per box.



SCANLAN® SOLEM® Vein Holder

Ideal for assisting surgeons during vessel anastomosis, the SCANLAN® Solem™ vein holder gently yet firmly grasps vein grafts and precisely positions suture without trauma. Its slim profile and concave tip assures optimal exposure and the lightly serrated clip provides steady, consistent pressure, minimizing graft movement.

- **Radiopaque**
- **Not made with natural rubber latex**

| Cat. No. | Description |
|----------|-------------|
|----------|-------------|

| | |
|----------|--------------------------------------|
| 1001-760 | Solem™ vein holder, 6 1/2" (16.5 cm) |
|----------|--------------------------------------|

2 per sterile package, 5 packages per box.

