

ATON™ | SURGICAL

FIBER OPTIC HEADLIGHT POWERED BY HATTERAS



SEE BETTER  
SEE FARTHER  
SEE WHY



SINCE  
1971

CELEBRATING  
FIVE DECADES  
OF ILLUMINATION

**AtoN™ SURGICAL**

**HIGH-INTENSITY SURGICAL FIBER OPTIC HEADLIGHT WITH ADJUSTABLE SPOT, POWER WITH BFW HIGH BRIGHT HATTERAS, S-MODEL HATTERAS OR CHROMAPRO**

- Illumination Ranges at 16" (40 cm):
  - > High Bright Hatteras ~220,000 LUX
  - > S-Model Hatteras ~185,000 LUX
  - > ChromaPRO ~70,000 LUX
- Adjustable Iris Spot Diameter: .5"-4" (1.27 - 10.16 cm) at a working distance of 16" (40 cm)
- Fiber Optic Cables available in 7' and 10' lengths
- Hygienic Optic: Easy to Clean Twist-Off Feature
- Compatibility: ACMI Port
- Optic Weight: 1.5 oz (43.4g)

**HATTERAS/AtoN SYSTEM**

- Hatteras LED High-Intensity Surgical Light Source
- AtoN High-Intensity Headlight
- Choice of: Fiji or Key West Headband
- 10' Cobra Coil Cable
- Mobile Stand

**AtoN HEADLIGHT SYSTEM PRODUCT CODES:**

- BFW 1200 AtoN Headlight with Fiji Headband and 10' Cobra Coil Fiber Optic Cable
- BFW 1252 AtoN Headlight with Key West Headband and 10' Cobra Coil Fiber Optic Cable
- BFW 1255 AtoN Headlight with Fiji Headband and 10' Standard Fiber Optic Cable
- BFW 1257 AtoN Headlight with Key West Headband and 10' Standard Fiber Optic Cable
- BFW 1260 AtoN Headlight with Fiji Headband and 7' Standard Fiber Optic Cable
- BFW 1262 AtoN Headlight with Key West Headband and 7' Standard Fiber Optic Cable
- BFW 9860 7' AtoN Standard Fiber Optic Cable
- BFW 9865 10' AtoN Standard Fiber Optic Cable
- BFW 9866 10' AtoN Cobra Coil Fiber Optic Cable



**EXCELLENT FOR USE IN THE FOLLOWING DEPARTMENTS**

**HOSPITAL**

- CARDIOSURGERY
- NEUROSURGERY
- ORTHOPEDIC SURGERY
- PEDIATRIC SURGERY
- SPINE SURGERY
- THORACIC SURGERY
- TRANSPLANT SURGERY
- VASCULAR SURGERY

**SURGERY CENTER**

**ORTHOPEDIC SURGERY**



[BFWinc.com](http://BFWinc.com)

800-717-4673 | +1-502-899-1808

1550 SAMCO ROAD

RAPID CITY, SD, 57702 USA

**ENGINEERED AND ASSEMBLED IN THE U.S.A.** © 2024 BFW, Inc.

ALL TECHNICAL SPECIFICATIONS ARE SUBJECT TO CHANGE WITHOUT NOTICE

**ATON TESTED AND CERTIFIED TO:** ANSI/AAMI ES60601-1 (2005 + C1:09 + A2:10), CAN/CSA C22.2 NO 60601-1 Edition 2 (2008), CE