

# CLAMP JAW COVERS



## SCANLAN® SOFT-GRIP™

### Fabric Clamp Covers

Special foam and fabric combination provides a soft, non-slip clamp surface.

- **Special material** – assures no fibrous residue in operating field
- **Versatile** – designed to fit a variety of instruments
- **Ease of application** – designed with convenient pull tabs
- **Adaptable** - may be trimmed to desired length
- **Radiopaque**
- **Not made with natural rubber latex**

Cat. No.	Description
1001-740	SOFT-GRIP™ fabric clamp covers Length 5" (127 mm), Width 3/8" (9.5 mm)

2 per sterile package, 10 packages per box.

## SCANLAN® SURG-I-PAW®

### Instrument Jaw Covers

Designed for Intestinal, Thoracic, Pediatric, and other specialty types of surgery.

- **Atraumatic** – converts standard clamps into atraumatic instruments
- **Deep serrations** – provides firm, atraumatic occlusion
- **Contoured design** – provides secure pressure - will not slip or twist
- **Radiopaque**
- **Not made with natural rubber latex**

Cat. No.	Description
1001-72	1/10" (2.5 mm), Red
1001-73	2/10" (5.0 mm), White
1001-74	3/10" (7.5 mm), Blue

2 per sterile package, 10 packages per box

