

XL Transport Containers

Transporting *da Vinci*® SP®, *Xi*® endoscopes, and *EndoWrist*® instruments appropriately and effectively is a difficult task. The InstruSafe® XL Transport Container and Lifting Tray are perfectly paired to soak and protect your instruments during transport from the O.R. to the SPD.



INSTRUSAFE



IN-PIC-20-D

XL TRANSPORT CONTAINER

XL Transport Container with Drain

Holds:

- *da Vinci SP* endoscopes and *EndoWrist* Instruments
- *da Vinci Xi* endoscopes and *EndoWrist* Instruments

34 L x 17 ¼ W x 6 H (in)

864 L x 438 W x 152 H (mm)

IN-PIC-20-D



IN-PIC-20

XL TRANSPORT CONTAINER

XL Transport Container

Holds:

- *da Vinci SP* endoscopes and *EndoWrist* Instruments
- *da Vinci Xi* endoscopes and *EndoWrist* Instruments

34 L x 17 ¼ W x 6 H (in)

864 L x 438 W x 152 H (mm)

IN-PIC-20

InstruSafe is a registered trademark of Summit Medical, an Innovia Medical Company.
da Vinci®, *EndoWrist*®, *Xi*® and *SP*® are all registered trademarks of Intuitive Surgical, Inc.

Lifting Tray

Transporting *da Vinci SP, Xi* endoscopes, and *EndoWrist* instruments appropriately and effectively is a difficult task. The InstruSafe XL Transport Container and Lifting Tray are perfectly paired to soak and protect your instruments during transport from the O.R. to the SPD.



INSTRUSAFE



IN-PIC-20-D-LT



IN-PIC-20-LT



Single Endoscope Configuration

Multiple Instrument Configuration

IN-PIC-34

XL TRANSPORT CONTAINER SYSTEM

XL Transport Container with Drain + Lifting Tray

Holds:

- Lifting Tray
- 1 *da Vinci Xi* or *SP* endoscope **OR** up to 8 *da Vinci Xi* or *SP EndoWrist* Instruments and accessories

Lifting Tray	Transport Container
27 ½ L x 13 ½ W x 3 ⅞ H (in)	34 L x 17 ¼ W x 6 H (in)
699 L x 343 W x 98 H (mm)	864 L x 438 W x 152 H (mm)

IN-PIC-20-D-LT

XL TRANSPORT CONTAINER SYSTEM

XL Transport Container + Lifting Tray

Holds:

- Lifting Tray
- 1 *da Vinci Xi* or *SP* endoscope **OR** up to 8 *da Vinci Xi* or *SP EndoWrist* Instruments and accessories

Lifting Tray	Transport Container
27 ½ L x 13 ½ W x 3 ⅞ H (in)	34 L x 17 ¼ W x 6 H (in)
699 L x 343 W x 98 H (mm)	864 L x 438 W x 152 H (mm)

IN-PIC-20-LT

LIFTING TRAY

Single Endoscope Configuration

Holds:

- 1 *da Vinci Xi* or *SP* endoscope and accessories

Multiple Instrument Configuration

Holds:

- Up to 8 *da Vinci Xi* or *SP EndoWrist* Instruments and accessories

27 ½ L x 13 ½ W x 3 ⅞ H (in)	699 L x 343 W x 98 H (mm)
------------------------------	---------------------------

IN-PIC-34

XL TRANSPORT CONTAINER

LID FEATURES

Built-in biohazard symbol and dedicated dry erase “soak start time” space on lid



*Top view of lid

SPILL REDUCTION

Inner gasket and strong silicone latches reduce spills

CART WASH SAFE

Highly durable material designed to withstand everyday use — cart wash safe — maximum temperature of 230°F

PROPER SOAKING ORIENTATION

Horizontally soak robotic instruments

VOLUME MEASURING

Liquid volume feature to accurately measure detergent solutions



QUICK DRAIN HOSE

Quick drain feature — no need to tip container to drain



*Features on top view of lid are emphasized for illustrative purposes only.

LIFTING TRAY

LOADING IDENTIFIERS

Permanent loading identifiers for correct and efficient instrument placement



LARGE HANDLES

Large handles on each end for safe and easy lifting

HIGHLY DURABLE

Withstands automated washing and use with enzymatic cleaners and detergents

ACCESSORY BASKET


Removable small accessory basket



CONFIGURATION EXAMPLES

Single Endoscope Configuration

Holds: 1 *da Vinci Xi* or *SP* endoscope



Multiple Instrument Configuration

Holds: up to 8 *da Vinci Xi* or *SP EndoWrist* Instruments

