

**NEW**

# Curved SCANLAN® RELIANCE Bulldog Clamp

The SCANLAN® Reliance Bulldog Clamp has been designed to provide secure, robust clamping pressure to assist in the temporary occlusion of vessels. Manufactured with the highest quality stainless steel, its 1 x 2 DeBakey jaws ensure a firm yet atraumatic clamping pressure. A unique suture port allows optional suture tagging to enhance surgical location and recovery.

## SCANLAN® Reliance Bulldog Clamps Curved

### Description

3795-53		<b>40 mm Curved DeBakey Jaws</b> 2 3/4" (7 cm)
---------	---	---

Firm Sale - May Be Available For Evaluation

**PLEASE NOTE:**

Not available in the European Union or Australia.

All products are not available in all regions. Please contact Scanlan customer service for product availability in your region.



## SCANLAN® Reliance Bulldog Applier / Introducer

		Shaft	Working Length	Total Length
9009-395	Straight Jaws	8 mm	13 3/4" (35 cm)	17 1/3" (44 cm)
9009-393	Straight Jaws	8 mm	5" (12.7 cm)	9" (23 cm)

### Features Include:

- 1 x 2 DeBakey jaws
- 350 grams closing pressure\*
- Suture port
- Highest quality stainless steel
- Reusable / durable
- Low profile design
- Radiopaque

### Reliance Bulldog Applier / Introducer

- Custom jaw allows for secure deployment of bulldog into the surgical field
- Designed for use with 12 mm cannula

\* Approximate clamping pressure

**For more information or to place an order,**  
please contact your local Scanlan representative, or call...

U.S. & Canada: 800-328-9458  
International: 651-298-0997  
FAX: 651-298-0018  
www.scanlaninternational.com



The World's Finest Surgical Instrumentation  
Since 1921