



Lithotomy Colon / Rectal



Abdominal Colon / Rectal



Prone Jack Knife Colon / Rectal



COLON / RECTAL SYSTEM #SL82002

REF	IMAGE	PART #	ITEM DESCRIPTION	QTY	APPLICATION
Tray 1					
A		43902AC	Elite III Rail Clamp with 1 Cam Joint 18"	1	Attaches to the rail of the OR table and holds the crossbar in position.
B		43900BC	24" Crossbar with 2 Cam III Joints (11" x 13")	1	Attaches to the rail clamp and holds the angled arms in place around the incision. May also support retractor handles.
C		44118N	18" Angled Arm (8" x 10") @ 45°	2	Attaches to the crossbar and holds the retractor handles in place around the incision.
		44124N	24" Angled Arm (8" x 16") @ 45°	2	
		44118B	18" Angled Closure Arm (9" x 9")	1	
		50000G	Instrument Case (23" x 11" x 3 1/2")	1	Holds instruments.
Tray 2					
D		SL42126WGP	Cam II Clip-on Quick Angle 8"	3	Attaches retractor blades to frame and angles retractor blades as needed, 45° up or down.
		SL45006CGP	Micro-Adjustable II Clip-on Quick Angle 10"	2	Attaches blades to frame and angles blades as needed, 45° up or down. Provides precise micro-adjustments.
		60020	T-Handle	2	May be used on Quick Angle Handle to angle blade 45° up or down.
		SL46143	Fenestrated 44mm x 89mm (1 3/4" x 3 1/2")	2	Abdominal or incisional wall retraction.
		SL46170	Deaver 25mm x 127mm (1" x 5")	2	Soft tissue retraction.
		SL46185	Hill-Ferguson 32mm x 76mm (1 1/4" x 3")	1	Exposure for hemorrhoid surgery.
		SL46186	Hill-Ferguson 25mm x 64mm (1" x 2 1/2")	1	
E		SL46187	Rectal Blade 13mm x 117mm (1/2" x 4 1/2")	2	Soft tissue retraction to expose bowel.
		SO46187A	St. Marks [no lip] 64mm x 178mm (2 1/2" x 7")	1	Soft tissue retraction.
		SO46188	St. Marks [with lip] 64mm x 178mm (2 1/2" x 7")	1	
		SO46189	St. Marks [with lip] 64mm x 140mm (2 1/2" x 5 1/2")	1	
		50000G	Instrument Case (23" x 11" x 3 1/2")	1	Holds instruments.