

SCANLAN®

Sundt™ Suction System

Rotatable & Non-Rotatable Graduated Suction

The SCANLAN® SUNDT™ Graduated Suction System provides complete control to atraumatically retract, probe, lift, and evacuate fluid or tissue from the surgical field. The SUNDT™ system offers a choice of seven tip diameters, each available in three lengths. The universal suction handle easily accommodates all tip diameters, lengths, and configurations.



Feel the Difference of Precision Craftsmanship



Highest Quality Surgical Instruments
Since 1921

SCANLAN®

Sundt™ Rotatable & Non-Rotatable Graduated Suction System

Features:

- Atraumatic tips may be used to retract or probe delicate tissue
- Precise fingertip control for atraumatic tissue manipulation or aspiration
- Each suction tip is tapered to reduce clogging
- One universal suction handle accommodates every tip configuration
- Tips and handle easily separate for cleaning and sterilization
- Seven suction tip diameters, ranging from 3 Fr. to 11 Fr.
- Three suction tip lengths
- Rotating or Non-rotating tips
- Guide pin securely holds tips in place at angle (for non-rotating tips)
- Universal SCANLAN® V.I.P.™ Sterilization tray holds 7 tips plus handles, fits all lengths
- Tips, handle, and tray may be purchased separately or in sets



SCANLAN®

Universal Suction Handle

Cat. No.	Description	Total Length
9009-930	Universal Suction Handle (Fits all tip lengths and french sizes)	3" (7.5 cm)



SCANLAN®

Suction Tips



Each Scanlan instrument is individually handcrafted to provide the surgeon with the highest quality, longest-lasting, precision instrument possible. Reflecting this unique process, all stated dimensions are approximate.

Not all products are available in all regions. Please consult your Scanlan Representative for products available in your region.

Please check with your Scanlan Representative for Firm Sale status. Custom modifications may be available upon request.



Sundt™ Rotatable & Non-Rotatable Graduated Suction System

SCANLAN®

Sundt™ Rotatable Sets

Each set has all seven suction tips and includes one universal suction handle and one V.I.P.™ sterilization tray.

Set Components	9009-940 Size: 4-1/2" (11.5 cm)	9009-960 Size: 5-1/4" (13.5 cm)	9009-980 Size: 7-7/8" (19.8 cm)
• 3 Fr. Tip	9009-941	9009-961	9009-981
• 4 Fr. Tip	9009-942	9009-962	9009-982
• 5 Fr. Tip	9009-943	9009-963	9009-983
• 6 Fr. Tip	9009-944	9009-964	9009-984
• 7 Fr. Tip	9009-945	9009-965	9009-985
• 9 Fr. Tip	9009-947	9009-967	9009-987
• 11 Fr. Tip	9009-949	9009-969	9009-989
Universal Suction Handle	9009-930	9009-930	9009-930
SCANLAN® V.I.P.™ Sterilization Tray	2136-03	2136-03	2136-03

SCANLAN®

Sundt™ Non-Rotatable Sets

Each set has all seven suction tips and includes one universal suction handle and one V.I.P.™ sterilization tray.

Set Components	9009-950 Size: 4-1/2" (11.5 cm)	9009-970 Size: 5-1/4" (13.5 cm)	9009-990 Size: 7-7/8" (19.8 cm)
• 3 Fr. Tip	9009-951	9009-971	9009-991
• 4 Fr. Tip	9009-952	9009-972	9009-992
• 5 Fr. Tip	9009-953	9009-973	9009-993
• 6 Fr. Tip	9009-954	9009-974	9009-994
• 7 Fr. Tip	9009-955	9009-975	9009-995
• 9 Fr. Tip	9009-957	9009-977	9009-997
• 11 Fr. Tip	9009-959	9009-979	9009-999
Universal Suction Handle	9009-930	9009-930	9009-930
SCANLAN® V.I.P.™ Sterilization Tray	2136-03	2136-03	2136-03

The SCANLAN® Sundt™ Graduated Suction System may be purchased as a set or as individual components.

Each Scanlan instrument is individually handcrafted to provide the surgeon with the highest quality, longest-lasting, precision instrument possible. Reflecting this unique process, all stated dimensions are approximate.

Not all products are available in all regions. Please consult your Scanlan Representative for products available in your region.

Please check with your Scanlan Representative for Firm Sale status. Custom modifications may be available upon request.



SCANLAN®

Sundt™ Rotatable & Non-Rotatable Graduated Suction System

SCANLAN®

V.I.P.™ Sterilization Tray



Cat. No.	Description	Total Length
2136-03	SCANLAN® V.I.P.™ Sterilization Tray, with Scan-Guard™ silicone brackets in base	6" x 9" x 2" (15 cm x 23 cm x 5 cm)

For more information or to place an order,
please contact your local Scanlan representative, or call...

U.S. & Canada: 800-328-9458 International: 651-298-0997
www.scanlaninternational.com

Each Scanlan instrument is individually handcrafted to provide the surgeon with the highest quality, longest-lasting, precision instrument possible. Reflecting this unique process, all stated dimensions are approximate.

Not all products are available in all regions. Please consult your Scanlan Representative for products available in your region.

Please check with your Scanlan Representative for Firm Sale status. Custom modifications may be available upon request.

