

BAYPORT WILL CHARGE UP
TO 6 BFW 3-CELL AND 4-CELL
LITHIUM-ION BATTERIES
(BATTERIES SOLD SEPARATELY)

—BFW 9020 BC6 —

SEE BETTER
SEE FARTHER
SEE WHY




SINCE 1971 CELEBRATING FIVE DECADES OF ILLUMINATION

BAYPORT™ 6-BAY BATTERY CHARGING SOLUTION

- On a Wall Mount Rail Kit, BayPort can simultaneously charge up to 6 batteries and store up to 4 headlights on a single module.
- On a Mobile Stand Kit, BayPort stores headlights and charges batteries in flexible combinations to suit. 4 modules required.
- On a Table Stand Kit, BayPort can conveniently store up to 2 headlights and charge up to 6 batteries.
- Lightning Fast charges - 80% of full power in 1 hour and 100% of full power in 3 hours.
- Exclusively designed to charge BFW 3- and 4-cell Lithium-ion batteries.
- Hook design will hold headlights, cable, and holster.
- Designed to store AtoN, Pharos, Dover, Bristol-Plus, and Daymark.
- Minor assembly required.
- Limited Manufacturer Warranty applies.

WALL MOUNT CONFIGURATION



—BFW 9020 BC6 + BFW 5440 MR
+ BFW 5440 TDK + BFW 5440 BDK—

TABLE STAND CONFIGURATION



—BFW 9020 BC6 + BFW 5440 TS—

MOBILE STAND CONFIGURATION



—BFW 9020 BC6 + BFW 5440 MS
+ BFW 5440 W2—

(MOBILE STANDS REQUIRE 4 MODULES.
HEADLIGHTS AND BATTERIES SHOWN
ARE SOLD SEPARATELY)

BAYPORT PRODUCT AND ACCESSORY CODES

BFW 9020 BC6	BayPort 6-bay Battery Charging Solution (Charges up to 6 BFW 3-cell and 4-cell Lithium-ion batteries exclusively)
BFW 5440 MR	BayPort Wall Mount Rail Kit
BFW 5440 MS	BayPort Mobile Stand Kit
BFW 5440 TS	BayPort Table Stand Kit
BFW 5440 TDK	BayPort Top Dual Hanger Kit
BFW 5440 BDK	BayPort Bottom Dual Hanger Kit
BFW 9893 PS	BayPort Mobile Stand Power Strip



BFWinc.com

800-717-4673 | +1-502-899-1808

1550 SAMCO ROAD

RAPID CITY, SD, 57702 USA

PLEASE NOTE

- Product specifications may be modified by the manufacturer at any time, descriptions for reference only.
- Configurations shown combine products and accessories that are sold separately. Consult your BFW specialist for the optimal configuration to suit your needs.

ENGINEERED AND ASSEMBLED IN THE U.S.A. © 2024 BFW, Inc.

ALL TECHNICAL SPECIFICATIONS ARE SUBJECT TO CHANGE WITHOUT NOTICE

BAYPORT 6-BAY CHARGING SOLUTION TESTED AND CERTIFIED TO: IEC 60601-1 EDITION 3, IEC 60601-1-2 Edition 4, ANSI/AAMI ES60601-1: A1:2012, C1:2009/(R)2012 and A2:2010/(R)2012, CSA CAN/CSA-C22.2 NO. 60601-1:14, IEC 60601-1-6:2010 (Third Edition) + A1:2013, CE