



SCANLAN®

# Microsurgical Instruments



**Feel the Difference of Precision Craftsmanship**

World's Finest Surgical Instruments  
Since 1921



# Super-Micro Fine Instrumentation

Highest Quality Materials and Manufacturing to Provide a Difference You Can Feel.

## Committed To Our Customers...

Dear Customer,

The Scanlan family has specialized in the design and manufacture of the world's finest surgical instrumentation and single-use products since 1921. Now in our fourth generation, our family continues its dedication to providing you with the best surgical products available in the world. This brochure is a collection of our precision products, all designed specifically to address your needs within the O.R. and S.P.D.

Our goal is to continuously enhance our designs to serve you, your hospital and ultimately, your patients. Thank you for your continued support and demand for the finest surgical products.

Sincerely,

Timothy M. Scanlan  
President and Chief Executive Officer

Russell T. Scanlan  
Chief Operating Officer



Joseph Scanlan

Russell Scanlan

Tim Scanlan

Adam Scanlan

Chenoa Scanlan

### SCANLAN® Super-Micro Fine

#### Needle Holders

- **High Quality Titanium and Stainless Steel**
- **Titanium Needle Holders are Non-magnetic** – eliminates needle jump
- **Titanium Needle Holders are Lightweight** – reduces the chance for wrist fatigue
- **Low Profile Tapered Shanks** for improved visibility under the microscope
- **Designed For Use With 9-0 Suture** (100 micron needles) or smaller
- **0.3 mm Super-Micro Fine tip**
- **DIAMOND DUST™** jaws provides a secure grip while protecting the needle from burring
- **Non-locking**



### SCANLAN® Super-Micro Fine

#### Forceps

- **High Quality Titanium and Stainless Steel**
- **Titanium Forceps are Non-magnetic** – eliminates needle from clinging to instrument
- **Titanium Forceps are Lightweight** – reduces the chance for wrist fatigue
- **Low Profile Tapered Shanks** for improved visibility under the microscope
- **Specifically designed for manipulation of small vessels**
- **0.3 mm Super-Micro Fine Tip**
- **DIAMOND DUST™** coating provides superior grasping of tissue



### SCANLAN® Super-Micro Fine

#### Scissors

- **High Quality Hardened Stainless Steel** for industry leading durability
- **SCANLAN® Premier™** scissors use proprietary technology that enhances the hardness and corrosion resistant properties of stainless steel
- **100% Cut Tested** to ensure a high level of quality
- **12-month Guarantee** on the sharpness of Scanlan® Premier™ scissors



All Scanlan Instruments are backed by a 10-year Guarantee against manufacturing material defects



# SCANLAN® MICROSURGICAL NEEDLE HOLDERS

## SCANLAN® Super-Micro Fine Needle Holders for 9-0 Suture



3303-646



3303-622

Cat. No.	Description	Metal	Handle	Recommended Suture Size	Total Length
<b>NEW</b> 3303-646	SCANLAN® Super-Micro Fine Titanium Needle Holder Curved, 0.3 mm DIAMOND DUST™ Jaws Regular Box Lock Without Ratchet	Titanium	Round	9-0 and Smaller	5" (12.7 cm)
3303-603	SCANLAN® Super-Micro Fine Titanium Needle Holder Straight, 0.3 mm DIAMOND DUST™ Jaws Regular Box Lock Without Ratchet	Titanium	Flat		
<b>NEW</b> 3303-664	SCANLAN® Super-Micro Fine Titanium Needle Holder Straight, 0.3 mm DIAMOND DUST™ Jaws, Regular Box Lock Without Ratchet	Titanium	Round		
<b>NEW</b> 3303-622	SCANLAN® Super-Micro Fine Titanium Needle Holder Straight, 0.3 mm DIAMOND DUST™ Jaws, Regular Box Lock Without Ratchet	Titanium	Flat		7-1/4" (18.5 cm)
3303-667	SCANLAN® Super-Micro Fine Titanium Needle Holder Curved, 0.3 mm DIAMOND DUST™ Jaws, Regular Box Lock Without Ratchet	Titanium	Round		

Each Scanlan instrument is individually handcrafted to provide the surgeon with the highest quality, longest-lasting, precision instrument possible. Reflecting this unique process, all stated dimensions are approximate.

Not all products are available in all regions. Please consult your Scanlan Representative for products available in your region.

Please check with your Scanlan Representative for Firm Sale status. Custom modifications may be available upon request.

For a complete list of models visit: [www.scanlaninternational.com](http://www.scanlaninternational.com)



# SCANLAN® MICROSURGICAL NEEDLE HOLDERS

## SCANLAN® Micro Needle Holders for 8-0 Suture



3003-177-4

Cat. No.	Description	Metal	Handle	Recommended Suture Size	Total Length
3003-177-3	SCANLAN® Titanium Jacobson Micro Needle Holder Curved, 0.3 mm DIAMOND DUST™ Jaws Regular Box Lock Without Ratchet	Titanium	Round	8-0 and Smaller	7-1/4" (18.5 cm)
3003-177-4	SCANLAN® Titanium Jacobson Micro Needle Holder Curved, 0.4 mm DIAMOND DUST™ Jaws Regular Box Lock Without Ratchet	Titanium			

## SCANLAN® Micro Needle Holders for 7-0 Suture



3003-207

Cat. No.	Description	Metal	Handle	Recommended Suture Size	Total Length
6006-44	SCANLAN® Castroviejo Micro Needle Holder Straight DIAMOND DUST™ Jaws Regular Box Lock With Ratchet	Stainless Steel	Flat	7-0 and Smaller	4-3/4" (12 cm)
3003-207	SCANLAN® Titanium Castroviejo Micro Needle Holder Straight DIAMOND DUST™ Jaws Regular Box Lock With Ratchet	Titanium			5-3/4" (14.5 cm)

Each Scanlan instrument is individually handcrafted to provide the surgeon with the highest quality, longest-lasting, precision instrument possible. Reflecting this unique process, all stated dimensions are approximate.

Not all products are available in all regions. Please consult your Scanlan Representative for products available in your region.

Please check with your Scanlan Representative for Firm Sale status. Custom modifications may be available upon request.

For a complete list of models visit: [www.scanlaninternational.com](http://www.scanlaninternational.com)




# SCANLAN® MICROSURGICAL NEEDLE HOLDERS

SCANLAN®

## Micro Needle Holders

for 6-0 Suture

3003-452



Cat. No.	Description	Metal	Handle	Recommended Suture Size	Total Length
3003-452	SCANLAN® Titanium Richards Micro Needle Holder Extra Fine Straight DIAMOND DUST™ Jaws With Ratchet	Titanium	Ring	6-0 and Smaller	6-3/4" (17 cm)

SCANLAN®

## Micro Needle Holders

for 5-0 Suture

4155-06



Cat. No.	Description	Metal	Handle	Recommended Suture Size	Total Length
4155-07	SCANLAN® Webster Needle Holder Straight delicate jaws Offset bows	Stainless Steel	Ring	5-0 and Smaller	5-1/8" (13 cm)
4155-06	SCANLAN® Webster Needle Holder Straight delicate jaws	Stainless Steel			5-1/8" (13 cm)
6006-04	SCANLAN® Dennis Delicate Needle Holder Straight jaws	Stainless Steel			7" (18 cm)
3003-422	SCANLAN® Microvascular Needle Holder Straight Jaws With Ratchet	Stainless Steel			7-1/8" (18.1 cm)

Each Scanlan instrument is individually handcrafted to provide the surgeon with the highest quality, longest-lasting, precision instrument possible. Reflecting this unique process, all stated dimensions are approximate.

Not all products are available in all regions. Please consult your Scanlan Representative for products available in your region.

Please check with your Scanlan Representative for Firm Sale status. Custom modifications may be available upon request.

For a complete list of models visit: [www.scanlaninternational.com](http://www.scanlaninternational.com)



# SCANLAN® MICROSURGICAL NEEDLE HOLDERS

SCANLAN®

## Micro Needle Holders

for 5-0 Suture

3003-90A



Cat. No.	Description	Metal	Handle	Recommended Suture Size	Total Length
3003-90A	SCANLAN® Titanium Castroviejo Micro Needle Holder Straight DIAMOND DUST™ Jaws Streamline Box Lock with Ratchet	Titanium	Round	5-0 and Smaller	4-3/4" (12 cm)
4155-06	SCANLAN® Webster Needle Holder Delicate DIAMOND DUST™ Jaws	Stainless steel	Ring		5-1/8" (13 cm)
6006-24	SCANLAN® Micro Needle Holder DIAMOND DUST™ Jaws	Stainless steel	Round		5-3/4" (14.5 cm)
6006-50	SCANLAN® Castroviejo Micro Needle Holder Straight DIAMOND DUST™ Jaws Regular Box Lock With Ratchet	Stainless steel			
6006-56-3	SCANLAN® Castroviejo Micro Needle Holder Curved, 0.3mm tapered DIAMOND DUST™ Jaws Regular Box Lock Without Ratchet	Stainless steel	Round		6" (15 cm)
3003-94-4	SCANLAN® Titanium Castroviejo Micro Needle Holder Curved, 0.4mm DIAMOND DUST™ Jaws Streamline Box Lock Without Ratchet	Titanium			
6006-28	SCANLAN® Micro Needle Holder DIAMOND DUST™ Jaws	Stainless steel	Ring	6-1/4" (16 cm)	

Each Scanlan instrument is individually handcrafted to provide the surgeon with the highest quality, longest-lasting, precision instrument possible. Reflecting this unique process, all stated dimensions are approximate.

Not all products are available in all regions. Please consult your Scanlan Representative for products available in your region.

Please check with your Scanlan Representative for Firm Sale status. Custom modifications may be available upon request.

For a complete list of models visit: [www.scanlaninternational.com](http://www.scanlaninternational.com)



# SCANLAN® MICROSURGICAL FORCEPS

## SCANLAN® Super-Micro Fine Forceps



3303-410



3303-470

Cat. No.	Description	Metal	Handle	Total Length
<b>NEW</b> 3303-452	SCANLAN® Super-Micro Fine Titanium Forceps Angled 0.3 mm DIAMOND DUST™ Tips	Titanium	Round	5" (12.7 cm)
<b>NEW</b> 3303-410	SCANLAN® Super-Micro Fine Titanium Forceps Straight 0.3 mm DIAMOND DUST™ Tips	Titanium	Flat	
3303-470	SCANLAN® Super-Micro Fine Titanium Forceps Straight 0.3 mm DIAMOND DUST™ Tips	Titanium	Round	7-1/4" (18.5 cm)
<b>NEW</b> 3303-472	SCANLAN® Super-Micro Fine Titanium Forceps Curved 0.3 mm DIAMOND DUST™ Tips	Titanium		
<b>NEW</b> 3303-430	SCANLAN® Super-Micro Fine Titanium Forceps Straight 0.3 mm DIAMOND DUST™ Tips	Titanium	Flat	7-1/4" (18.5 cm)
3303-434	SCANLAN® Super-Micro Fine Titanium Forceps Angled 0.3 mm DIAMOND DUST™ Tips	Titanium	Flat	

Each Scanlan instrument is individually handcrafted to provide the surgeon with the highest quality, longest-lasting, precision instrument possible. Reflecting this unique process, all stated dimensions are approximate.

Not all products are available in all regions. Please consult your Scanlan Representative for products available in your region.

Please check with your Scanlan Representative for Firm Sale status. Custom modifications may be available upon request.

For a complete list of models visit: [www.scanlaninternational.com](http://www.scanlaninternational.com)



# SCANLAN® MICROSURGICAL FORCEPS

## SCANLAN® Micro Forceps



3003-154-3

Cat. No.	Description	Metal	Handle	Total Length
4004-03-3A	SCANLAN® Micro Suture Tying Forceps 45 Degree Angled 0.3 mm, DIAMOND DUST™ Tips	Stainless Steel	Round	5-1/2" (14 cm)
3003-154-3	SCANLAN® Micro Suture Tying Forceps Gerald Type Taper Straight 0.3 mm DIAMOND DUST™ Tips	Titanium		6" (15 cm)
3003-156-3	SCANLAN® Micro Suture Tying Forceps Gerald Type Taper Curved 0.3 mm DIAMOND DUST™ Tips	Titanium		6" (15 cm)
3003-158-3	SCANLAN® Micro Suture Tying Forceps Gerald Type Taper Straight 0.3 mm DIAMOND DUST™ Tips	Titanium		7-1/4" (18.5 cm)
3003-159-3	SCANLAN® Micro Suture Tying Forceps Gerald Type Taper Curved 0.3 mm DIAMOND DUST™ Tips	Titanium		7-1/4" (18.5 cm)
4004-97-3A	SCANLAN® Micro Suture Tying Forceps 45 Degree Angled 0.3mm, DIAMOND DUST™ Tips	Stainless Steel		7-1/4" (18.5 cm)

## SCANLAN® Specialty Forceps



4004-29

Cat. No.	Description	Metal	Handle	Total Length
4004-05	SCANLAN® Jewelers Forceps Straight Sharp Smooth Tips	Stainless Steel	Flat	4-1/4" (11 cm)
4004-07	SCANLAN® Jewelers Forceps Straight Sharp DIAMOND DUST™ Tips	Stainless Steel		
4004-29	SCANLAN® Adson Forceps 1.2 mm Tips, 1x2 Teeth	Stainless Steel	Flat	4-3/4" (12 cm)
4004-19	SCANLAN® Jewelers Forceps Straight sharp DIAMOND DUST™ Tips	Stainless Steel		5-1/8" (13 cm)
4004-63	SCANLAN® McIndoe Forceps Serrated Tips, Fine 1x2 Teeth	Stainless Steel		6-1/4" (16.5 cm)

Each Scanlan instrument is individually handcrafted to provide the surgeon with the highest quality, longest-lasting, precision instrument possible. Reflecting this unique process, all stated dimensions are approximate.

Not all products are available in all regions. Please consult your Scanlan Representative for products available in your region.

Please check with your Scanlan Representative for Firm Sale status. Custom modifications may be available upon request.

For a complete list of models visit: [www.scanlaninternational.com](http://www.scanlaninternational.com)



# SCANLAN® MICROSURGICAL SCISSORS

SCANLAN®

## Super-Micro Fine Spring Style Scissors



7007-712



7007-768

Cat. No.	Description	Metal	Handle	Total Length
<b>NEW</b> 7007-712	SCANLAN® <i>Premier</i> ™ Super-Micro Fine Jacobson Scissors Curved Nano Blades	Stainless Steel	Flat	5" (12.7 cm)
<b>NEW</b> 7007-729	SCANLAN® <i>Premier</i> ™ Super-Micro Fine Jacobson Scissors Curved Nano Blades	Stainless Steel	Round	7-1/4" (18.5 cm)
<b>NEW</b> 7007-768	SCANLAN® <i>Premier</i> ™ Super-Micro Fine Jacobson Scissors Straight Nano Blades	Stainless Steel		
7007-124	SCANLAN® Jacobson Micro Scissors Curved Fine, Sharp Tips	Stainless Steel	Flat	7" (18 cm)
7007-120	SCANLAN® Jacobson Micro Scissors Straight Fine, Sharp Tips	Stainless Steel		
7007-140	SCANLAN® Jacobson Vannas Micro Scissors 45° Angled on Flat Micro Tips	Stainless Steel		

Each Scanlan instrument is individually handcrafted to provide the surgeon with the highest quality, longest-lasting, precision instrument possible. Reflecting this unique process, all stated dimensions are approximate.

Not all products are available in all regions. Please consult your Scanlan Representative for products available in your region.

Please check with your Scanlan Representative for Firm Sale status. Custom modifications may be available upon request.

For a complete list of models visit: [www.scanlaninternational.com](http://www.scanlaninternational.com)



# SCANLAN® MICROSURGICAL SCISSORS

SCANLAN®

## Micro Spring Style Scissors



7007-128

Cat. No.	Description	Metal	Handle	Total Length
7007-404	SCANLAN® <i>Premier</i> ™ Castroviejo Micro Scissors 45° Angled Blades, 1 Blade with Bead	Stainless Steel	Flat	4-1/8" (10.5 cm)
7007-108	SCANLAN® Castroviejo Micro Scissors Curved Blades, Sharp Tips	Stainless Steel	Round	5-1/2" (14 cm)
7007-538	SCANLAN® <i>Premier</i> ™ Jacobson Micro Scissors Straight Sharp, Nano Blades	Stainless Steel	Flat	6" (15 cm)
7007-481	SCANLAN® Jacobson Micro Scissors 45° Angled Blades, 1 Blade with Bead tip	Stainless Steel	Flat	6-1/2" (16.5 cm)
7007-128	SCANLAN® Jacobson Micro Scissors Straight Blades, Sharp Tips	Stainless Steel	Round	7" (18 cm)
7007-130	SCANLAN® Jacobson Micro Scissors Curved Blades, Sharp Tips	Stainless Steel	Round	7" (18 cm)

Each Scanlan instrument is individually handcrafted to provide the surgeon with the highest quality, longest-lasting, precision instrument possible. Reflecting this unique process, all stated dimensions are approximate.

Not all products are available in all regions. Please consult your Scanlan Representative for products available in your region.

Please check with your Scanlan Representative for Firm Sale status. Custom modifications may be available upon request.

For a complete list of models visit: [www.scanlaninternational.com](http://www.scanlaninternational.com)



SCANLAN®  
MICROSURGICAL SCISSORS

SCANLAN®  
Micro Ring Handle Scissors



7007-306SC

Cat. No.	Description	Metal	Handle	Total Length
7007-282SC	SCANLAN® Iris-Eye Scissors Straight Super Cut™ Blades, Sharp Tips, Black Platinum Ring Handles	Stainless Steel	Ring	4" (10 cm)
7007-210SC	SCANLAN® Stevens Tenotomy Scissors Curved Super Cut™ Blades, Black Platinum Ring Handles	Stainless Steel		4-5/8" (11.8 cm)
7007-214SC	SCANLAN® Baby Metzenbaum Scissors Curved Super Cut™ Blades, Black Platinum Ring Handles	Stainless Steel		4-3/4" (12 cm)
7007-272SC	SCANLAN® Ragnell Dissecting Scissors Curved Super Cut™ Blades, Flat Tips, Black Platinum Ring Handle	Stainless Steel		6" (15 cm)
7007-292SC	SCANLAN® C.V.H. Scissors Straight Super Cut™ Blades, Sharp Tips, Black Platinum Ring Handles	Stainless Steel		6" (15 cm)
7007-211SC	SCANLAN® Stevens Tenotomy Scissors Curved Super Cut™ Blades Black Platinum Ring Handles	Stainless Steel		7" (18 cm)
7007-218-1SC	SCANLAN® Metzenbaum Scissors Curved Super Cut™ Blades, Fine Tips, Black Platinum Ring Handles	Stainless Steel		
7007-306SC	SCANLAN® Face-Lift Scissors Curved Super Cut™ Blades, Semi-Sharp Outer Blade, Black Platinum Ring Handles	Stainless Steel		7-3/4" (19.5 cm)

Each Scanlan instrument is individually handcrafted to provide the surgeon with the highest quality, longest-lasting, precision instrument possible. Reflecting this unique process, all stated dimensions are approximate.

Not all products are available in all regions. Please consult your Scanlan Representative for products available in your region.

Please check with your Scanlan Representative for Firm Sale status. Custom modifications may be available upon request.

For a complete list of models visit: [www.scanlaninternational.com](http://www.scanlaninternational.com)



SCANLAN®  
CLAMPS / HAEMOSTATS

SCANLAN®  
Clamps / Haemostats



5005-64

Cat. No.	Description	Metal	Total Length
9009-148	SCANLAN® Ochsner Forceps Ring Handle, 90° Angle, Very Fine Jaws	Stainless Steel	4-1/2" (11.5 cm)
5005-02	SCANLAN® Vascular (Mosquito) Forceps Ring Handle, Delicate Curved Jaws	Stainless Steel	4-3/4" (12 cm)
5005-64	SCANLAN® Bailey Forceps Ring Handle, 90° Angled Very Fine Jaws	Stainless Steel	7" (18 cm)

Each Scanlan instrument is individually handcrafted to provide the surgeon with the highest quality, longest-lasting, precision instrument possible. Reflecting this unique process, all stated dimensions are approximate.

Not all products are available in all regions. Please consult your Scanlan Representative for products available in your region.

Please check with your Scanlan Representative for Firm Sale status. Custom modifications may be available upon request.

For a complete list of models visit: [www.scanlaninternational.com](http://www.scanlaninternational.com)



# SCANLAN® DILATORS

SCANLAN®

## Selber Graduated Dilator

Developed in cooperation with Jesse C. Selber, M.D., M.P.H.



Enlarged to Show Detail

2950-01  
Actual Suze

SCANLAN®

## Graduated Dilator



Curved dilator allows for placement in tighter anatomy

2950-02  
Actual Suze

Enlarged to Show Detail

Cat. No.	Description	Metal	Total Length
2950-01	SCANLAN® Selber Graduated Dilator Straight	Stainless Steel	5-1/2" (14 cm)
2950-02	SCANLAN® Graduated Dilator Curved	Stainless Steel	

SCANLAN®

## Micro Vessel Dilators



4004-150



4004-155

Cat. No.	Description	Metal	Total Length
4004-150	SCANLAN® Micro Vessel Dilator Spring Style, Flat Handle, Straight Tips	Stainless Steel	4-1/4" (11 cm)
4004-155	SCANLAN® Micro Vessel Dilator Spring Style, Round Handle, Straight Tips	Stainless Steel	4-7/8" (12.4 cm)

Each Scanlan instrument is individually handcrafted to provide the surgeon with the highest quality, longest-lasting, precision instrument possible. Reflecting this unique process, all stated dimensions are approximate.

Not all products are available in all regions. Please consult your Scanlan Representative for products available in your region.

Please check with your Scanlan Representative for Firm Sale status. Custom modifications may be available upon request.

For a complete list of models visit: [www.scanlaninternational.com](http://www.scanlaninternational.com)



# SCANLAN® SINGLE-USE BULLDOG CLAMPS

SCANLAN®

## VASCU-STATT® Bulldog Appliers



Cat. No.	Description	Metal	Total Length
1001-506	SCANLAN® VASCU-STATT® Applier Straight	Stainless Steel	8-1/2" (21.5 cm)
1001-506A	SCANLAN® VASCU-STATT® Applier Angled		
1001-506M	SCANLAN® Micro VASCU-STATT® Applier Straight		

SCANLAN®

## VASCU-STATT® Bulldog Clamps



Cat. No.	Description	Clamping Pressure
1001-507	SCANLAN® VASCU-STATT® Approximator Mini, Straight (Red)	Regular
1001-508	SCANLAN® VASCU-STATT® Approximator Mini, Angled 45° (Red)	Regular
1001-540	SCANLAN® VASCU-STATT® II Approximator Straight (Green)	Firm
1001-499	SCANLAN® Micro VASCU-STATT® Micro, Straight (Purple)	Regular
1001-529	SCANLAN® Micro VASCU-STATT® II Micro, Straight (Black)	Firm

Each Scanlan instrument is individually handcrafted to provide the surgeon with the highest quality, longest-lasting, precision instrument possible. Reflecting this unique process, all stated dimensions are approximate.

Not all products are available in all regions. Please consult your Scanlan Representative for products available in your region.

Please check with your Scanlan Representative for Firm Sale status. Custom modifications may be available upon request.

For a complete list of models visit: [www.scanlaninternational.com](http://www.scanlaninternational.com)



SCANLAN®

# MICROSURGICAL Instrumentation

Performance, quality and exceptional tactile response during the most complex procedures.

For more information or to place an order,  
please contact your local Scanlan representative,  
or call...

U.S. & Canada: 800-328-9458

International: 651-298-0997

FAX: 651-298-0018

[www.scanlaninternational.com](http://www.scanlaninternational.com)

