

# TIP-GUARD™

## Instrument Protectors

SCANLAN® TIP-GUARD™ - Ensures safe, cost-effective protection of blade tips and edges for all surgical instrumentation.

- **Economical** – provides long term use and is cost-effective
- **Superior protection** – safeguards delicate tips and edges of surgical instrumentation
- **Tested safe** – clinically tested for steam, flash and EO sterilization (sterility report available upon request)
- **Versatile** – conforms to various shapes and angles
- **Design variety** – available in colored, tinted, clear and vented
- **Radiopaque**
- **Not made with natural rubber latex**



# TIP-GUARD™

## Instrument Protectors

Colored



1001-89  
Black



1001-88  
Brown



1001-87  
Orange



1001-91  
Maroon



1001-90  
Yellow



1001-85  
Red



1001-86  
Green



1001-93  
Blue



1001-92  
White

1001-84 Variety Pack - 9 assorted color-coded sizes, 100 pieces per package

Tinted



1001-305  
Black



1001-304  
Brown



1001-303  
Orange



1001-307  
Maroon



1001-306  
Yellow



1001-301  
Red



1001-302  
Green



1001-309  
Blue



1001-308  
White

1001-300 Variety Pack - 9 assorted color-coded sizes, 100 pieces per package

Clear



1001-205



1001-204



1001-203



1001-207



1001-206



1001-201



1001-202



1001-209



1001-208

1001-200 Variety Pack - 9 assorted color-coded sizes, 100 pieces per package

Vented



1001-405  
Black



1001-404  
Brown



1001-403  
Orange



1001-407  
Maroon



1001-406  
Yellow



1001-401  
Red



1001-402  
Green



1001-409  
Blue



1001-408  
White

1001-400 Variety Pack - 9 assorted color-coded sizes, 100 pieces per package



<b>Width:</b>	1/8" (3.2 mm)	1/16" (1.6 mm)	1/16" (1.6 mm)	13/32" (10.3 mm)	3/16" (4.8 mm)	1/8" (3.2 mm)	7/64" (2.8 mm)	5/64" (2 mm)	1/16" (1.6 mm)
<b>Length:</b>	1" (25.4 mm)	1" (25.4 mm)	1" (25.4 mm)	3/4" (19 mm)	1" (25.4 mm)	1" (25.4 mm)	3/4" (19 mm)	3/4" (19 mm)	3/4" (19 mm)
<b>Height:</b>	1" (25.4 mm)	5/8" (15.9 mm)	3/8" (9.5 mm)						

